


The **SICrom[™]** Accelerated Liquid Chromatograph System **Onyx[™]** Advantages & Benefits: Flow Rates and Low Backpressure

The very high porosity of Onyx[™] allows you to run the column at very high flow rates. The figure below illustrates how the backpressures of both particle-based columns and monolithic columns are affected by flow rates.

Backpressure Comparison vs. Flowrates



Note: FIAlab Instruments is not affiliated with [Phenom](#)
[Return to the SICrom](#)

Please email or phone FIAlab Instruments for additional product information.
Email: fialab@flowinjection.com, Voice: 425-376-0450 or 1-800-963-1101, Fax: 425-376-0451
