

Sequential Injection for...

FIAlab

Leaders in Flow Injection Technology

On-line Process Monitoring

Sequential Injection (SI) is a liquid handling technique that can automate important analytical processes like **sample pretreatment** and **sample measurement**.

Applications include:

- **Nutrient Monitoring**
- **Product Characterization**
- **Industrial Mixing Processes**
- **Contaminant Monitoring**
- **and more...**

Operated by

SIASOFT

FOR WINDOWS 7, 8, 10

The FIAlab Process Analyzer is a revolutionary solution to for many analytical obstacles. It's versatility distinguishes it from any other biopharmaceutical tool.

Functions Include:

Dilution
Solid Phase Extraction
Heating
Gas Diffusion
Dialysis



FIAlab Process Analyzer