


# SIAlock



## SIAlock sterile barrier

ASK FOR A QUOTE



## ASEPTIC SAMPLING PRIOR TO LC-MS PRETREATMENT

Maintain sterile barrier with FIA Lab's SIAlock and SIAlock hub. With a maximum of 8 bioreactors per hub and 2 hubs per system.



Internal Volume closed	Internal Volume open	Dead Volume (uL)	Lifetime (# of cycles)
85 uL	185 uL	0 uL	10 million

## SIALOCK AND HUB SPECIFICATIONS

	Height	Width	Depth	Weight	Voltage	Current
<b>SIAlock</b>	8 cm	16 cm	23 cm	2 kg	24V	2.5A
<b>Hub</b>	31 cm	21 cm	29 cm	5 kg	9V-12V	2.5A